

Contact Lens Centre Product Guide



Contents

Custom Lenses
Eyecare Range
Solutions Range
Accessories Range



N X V

Custom Lenses

M G T H



F O E

• • •

RGP Material List

Material	Colours	WA	DK	UV	Non – UV
Acuity 18	LB B GN GY		21	✓	
Acuity 58	LB B GN GY		58	✓	
Acuity 100	B GN IB V C	23	111	✓	
Acuity 200 (small diam)	GN B C		200	✓	
Acuity 200 (large diam)	C B		200	✓	
Alberta 2	B GN GY C B	22	21		✓
Aquasil	B GN	35	50		✓
Bio GP 100	B GN	53	100		✓
Bio GP 100 (large diam)	GN	53	100		✓
Boston EO	B GY GN BN EB	48	82	✓	
Boston Equalens	B	30	72	✓	
Boston ES	B GY GN			✓	
Boston XO	B GN V IB	49	100	✓	
Boston XO (large diam)	C B	49	100	✓	
Boston XO2	GY	38	141	✓	
Fluoperm 30	GN C	13	30		✓
Fluoperm 60	UVGN	15	60		✓
Hybrid FS	B C GY		23		✓
Hybrid FS Plus	B GY GN		45		✓
Menicon Z*	LB		163	✓	
Onsi 56	B GN GY DB GN2 C V DKBN		56	✓	
Optimum Comfort	B GY GN LB	6	65	✓	
Optimum Extra	B GY GN LB	3	100	✓	
Optimum Extreme	B GY GN LB	6	125	✓	
Quantum II	B	49	100		✓

*Extra Charge

Colour **AM-** Amber **B-** Blue **BN-** Brown **C-** Clear **CB-** Crystal Blue
DB – Dark Blue **EB –**Electric Blue **GN –**Green **DKBN –** Dark Brown
GN2 – Forrest Green **GY –**Grey **GY2 –** Dark Grey **IB –** Ice Blue
LB – Light Blue **R –** Red **V –** Violet **Y-** Yellow

WA: Wetting Angle: the lower the angle, the more wettable the material

DK: The measurement of materials ability to pass oxygen; the higher the Number the more oxygen the materials can pass.

Soft Lens Materials

Materials	DK	Water Content %
Contaflex 38	9	38
Contaflex 58	21	58
Contaflex 67	30	67
Contaflex 75	43	75
Contaflex 77	39	77
Definitive 74	60	74
GM3	16	49
G M Advance 49	16	49
G M Advance 58	26	58
Benz G5X	22	59

DESIGNS

Bicurve, Aspheric, Kerasoft,
Centra ST, Q Toric, G Toric, Soft K, Rose K Soft,
Cantor & Nissel

Hybrid Lenses

- RSS



Eyecare Range



TheraTears Drops

- For moderate to severe dry eyes
- Contains 25 x 0.6ml units
- Active Ingredient: Sodium Carboxymethylcellulose 0.25%
- Preservative-free
- Lowest osmolarity (salt concentration) of any other Eye drops (181 mOsm/L)
- Recappable
- Patented electrolyte balance formula
- Great with soft and RGP contact lenses



TheraTears Liquid Gel

- For moderate to severe dry eyes
- Contains 30 x 0.6ml units
- Active ingredient: Sodium Carboxymethylcellulose 1.0%
- Preservative-free
- Lowest osmolarity (salt concentration) of any other Eye drops (181 mOsm/L)
- Recappable
- Patented electrolyte balance formula
- Can use during day or night-time
- Great with RGP and Ortho-K Lenses



TheraTears Sterilid Eyelid Cleanser

- For blepharitis and general eyelid hygiene
- Anti-bacterial
- Contains linalool oil (extract of tea tree oil)
- 48ml bottle
- Pumps out foam
- No expiry dates
- Can be used to help control bacteria on eyelids
With conditions like Conjunctivitis
- Great for pre-surgery to reduce the risk of eye infections
- Use in operating room if patient is allergic to iodine
- Can remove makeup



TheraTears Sterilid Antimicrobial Eyelid Cleanser

- Antimicrobial cleanser kills 99.9% of bacteria
- Contains Hypochlorous acid 0.01%
- Professional formula for superior daily eyelid and eyelash Cleansing.
- Effectively and quickly eliminates microorganisms from the eyelids
- pH balanced formula that is gentle on eyelids
- No rinse-soap free



Thera Tears Nutrition

- Dry eye supplement
- Long term treatment of dry eye, meibomian gland dysfunction (MGD) and blepharitis
- Helps naturally produce more tears and oils to the tear film
- Fish oil, lignan-free flaxseed oil & vitamin E
- 90 capsules per bottle (take 3 per day)
- For best results continue beyond 4 to 8 weeks (then maybe able to drop to 2capsules per day)
- Clinically proven



Evolve Hypromellose 0.3%

- For mild to moderate dry eye
 - 10ml bottle
 - 250 drops
 - 3-month use-up period after opening
 - Prevents damage to the cornea
- With patients who suffer from Keratoconjunctivitis sicca
- Preservative-free
 - Blue tip technology for accuracy
 - Squeezable bottle
 - Can be used with contact lenses



Please contact our office for further information



Evolve Carmellose 0.5%

- For mild to severe dry eye
- 10ml bottle
- 250 Drops
- 3-month use-up period after opening
- Osmoprotective
- Preservative-free
- Blue tip technology for accuracy
- Squeezable bottle
- Can be used with contact lenses



Please contact our office for further information



Evolve Eyelid Wipes

- Daily ocular cleansing and hygiene
- 20 Wipes
- Suitable for babies, sensitive skin, and adults
- Preservative-free
- Detergent Free



Please contact our office for further information



Oculocin Propo

- Relieves tired and irritated eyes
- Washes and hydrates eyes
- 10x0.5ml vials or 10ml bottle
- Active ingredients: Aloe Vera, Propolis and Chamomilla
- Preservative Free



Solutions Range



Oté Clean

- Optimal daily cleaning solution for hard and soft contact lenses
- Alcohol Based
- Ideal for plasma coated lenses
- Suitable for all kinds of RGP and soft contact lens materials Including SiHy.
- 40ml bottle



Eyeeye Saline

- PH: 7
- Low amount of preservatives
- Ingredients: Isotonic Sodium Chloride (Natriumchlorid 0.9%), EDTA 0.1%
- 360ml bottle
- Great for rinsing RGP or soft contact lenses
- Once opened discard after 1 month



Oté Clean/ Eyeeye Saline Combo

- 1x Oté Clean and 1x Eyeeye Saline



Menicon Menicare Plus

- RGP multipurpose solution
- Ingredients: Polyhexamethylene Biguanide 0.0005%, Poloxamer 0.5%, Hypromellose 0.275%, EDTA, buffered excipient
- 250ml bottle & 50ml bottle
- 1 free barrel lens case
- Effective against Acanthamoeba (cysts and trophozoites) in 2 hours
- Once opened discard after 3 months



Menicon Spray and Clean

- 15ml bottle
- Spray, rub and rinse GP lenses after removal at the end of every day to ensure lenses stay clean and comfortable from day one.
- Preservative free spray that is non-abrasive and compatible with all GP lenses

Designed to be used in conjunction with MeniCare Plus multipurpose solution and progent intensive weekly cleaner



Menicare Plus Combo

- 3 x Menicare Plus



Menicon Value Pack

- 3 x Menicare Plus, 1x Progent, 3x Barrel Lens Cases



Menicon Progent

- RGP protein removal solution
- 7 treatments (no case provided)
- Ingredients: A vial- Sodium Hypochlorite, B vial- Potassium Bromide
- Soak RGP lenses for 30 minutes then remove and rinse well (soaking for longer may discolour lenses)
- Effective in 5 minutes against bacteria, fungi, virus, Acanthamoeba and against prion in 30 minutes
- Can use weekly or as required



Menicon Progent Case

- Used for soaking RGP lenses in Progent or Menilab




Menicon Menilab

- In house RGP soaking solution
- Cleans, removes deposits, protein and sterilises RGP lenses
- Ingredients: Active chlorine solution qs 0.5% or 5000ppm of active chlorine, purified water qs 100%
- 250ml bottle
- Only 5 minutes soaking to sterilise lens to reduce risk on contamination of RGP trial lenses
- Soak RGP lenses for up to 30 minutes
- Always rinse well with Menicare Plus after soaking
- Once opened discard after 3 months



Accessories Range



DMV Range

- Come in packs of 10 (can order any quantity)
- DMV Classic (flat with a hole)
- DMV Ultra (flat with no hole)
- DMV 45 (angled without hole)
- DMV Magic Touch (angled with hole)
- DMV SLH (soft lens handler)
- DMV Scleral Cups (large flat with hole)
- DMV Vented Scleral Cups (large flat with vent all the way through, for inserting hybrid lenses)
- For removal and inserting contact lenses

